



- PERFECT FOR MOISTURE ANALYSIS
- PREVENTS SPLASHING OF LIQUID SAMPLES
- FACILITATES HEAT DISTRIBUTION
- BINDER LESS GLASS FIBER MATERIAL
- MEETS AOAC CRITERIA
- MANUFACTURED IN THE USA
- PAN SIZE: 3.5"

PART NUMBER	CONTAINS	DIMENSIONS
11334	50 / Box Pure Sample Aluminum Pans	90 mm ID / 110 mm OD
11336	200 / Box Pure Sample Glass Fiber Pads	90 mm OD / 0.015" thick